



CCD Based Multi-Wavelength X-ray Diffractometer

Keith Gendreau (GSFC), Zaven Arzoumanian (GSFC), Vanderlei Martins (GSFC), Rosemary Millham (OSU, SP systems), George Ricker (MIT), Peter Ford (MIT), Richard Starr (GSFC), Jack Trombka (GSFC)

Looking to expand




The Concept (for NASA Missions)

- A combination X-ray Diffractometer/Fluorescence instrument that requires **NO sample preparation**
- Definitively identifies minerals, **including water ice.**
- Provides particle/grain sizes, shapes, and orientations -> **clues to formation history**
- Based on flight proven X-ray CCD technology
- Has option to study aerosol mineral composition



Implementation






Application to VSE and Lunar Exploration



Water Ice Diffraction Features






Not the same as CHEMIN



GODDARD SPACE FLIGHT CENTER





$$\lambda = 12400 \text{ eV}\text{\AA}/E(\text{eV})$$



$$\theta = \text{Geometric function}(X)$$

$$n\lambda = 2d \sin(\theta)$$

**CXRDF
GSFC
Approach**



**Single
Wavelength
XRDF
Like CHEMIN**





Our new concept expands X-ray diffraction into a new dimension.



Prototype Data (Al 6061)




In a plot of energy vs 2Theta we see fluorescence lines as horizontally stretched events. The arcs in this space are due to diffraction and each arc represents a single d-spacing.



We can convert the energy of each event into a wavelength, take into account the geometry, and convert 2Theta and energy into a d-spacing value for each event using Bragg's Law. Now we can plot an orthogonal space of Energy vs. d-spacing.

The plot directly shows you XRF and XRD information.





Prototype Data




GODDARD SPACE FLIGHT CENTER

**Random
Piece of
Quartz
Found
Around
Building 2
at GSFC**


**Clearly see Si
Fluorescence
And ~10 XRD
Signatures.**

**Note that the
arcs are broken
up.**






Prototype Data (Quartz Sample)




2θ

The energy (wavelength) of the events at constant x are the same. As x gets smaller, the energy (wavelength) of the events gets higher (shorter) in accordance with Bragg's Law.



y



The combination of physical collecting area and broad wavelength coverage with spectroscopy makes the CXRDF insensitive to sample preparation.




Single Crystal of Calcium Carbonate


Here, we took a single crystal of calcium carbonate and oriented it 5 different ways by tipping and tilting at the sample position.

The image to the lower right shows how the diffracted spot moves around the focal plane. As the image is deflected to the left in the image, the energy of the photons diffracted shifts upward.

The image to the upper right shows that all these points line up on a “constant d-spacing arc” in the energy vs X plane.

A traditional XRD using a single wavelength could not do this.








Prototype Data



Hematite



Transform event list






Prototype Data



Water Ice:

*1st light
water ice
with 2nd
prototype
system*

*~ 1 minute of
data with ice
cube pressed
up to window*




*Beryllium
Instrument
Feature*

*Expected
water ice
d-spacing
arcs*



GODDARD SPACE
FLIGHT CENTER




27 June 2006 Keith Gendreau (GSFC)

ESTC 2006



GODDARD SPA
CE FLIGHT CEN
TER




27 June 2006 Keith Gendreau (GSFC)

ESTC 2006



GODDARD SPACE
FLIGHT CENTER




27 June 2006 Keith Gendreau (GSFC)

ESTC 2006




**Sensor Head Assembly:
~<5 cm on a side. Shielded**


X-ray Source



Small aperture (Possibly Be window)
For placing samples against



**Electronics: CCD Controller, Event Extractor,
Computer Interface or Event Logger**





Summary and Status



- A material analysis tool that requires no sample preparation and has potential to identify many minerals including water ice.
- Concept is simple
- No moving parts
- Prototype instruments made from commercial components
- Our MIT collaborators have been flying X-ray CCDs for 13+ years
- Estimated Mass/Power for new flight instrument using latest CCD electronic designs: <2 kg ; <2 watts average





Backup






XRD of Crystalline Structures



$$n \lambda = 2d \sin (\theta)$$




Aerosol Sampling Option



Optical probes for scattering and absorption (provides complementary information on size, refractive indices, etc.)



Instrument Description





X-ray Cameras from MIT CCD Group


MIT Center for Space Research has developed and fabricated X-ray cameras and flight electronics optimized for MIT/LL CCDs flown on the ASCA, Chandra, Astro-E, and HETE-2 missions.



ASCA/SIS (1993)




Chandra/ACIS (1999)



Astro-E/XIS (2000)




HETE SXC (2000)



CSR






Application to VSE and Lunar Exploration



- Definitive Identification of Minerals
 - Definitive Identification of Water-Ice
- Could Provide Insight into Lunar History via Analysis of structure of the d-spacings and particle/grain size, shape, and orientation distributions



GODDARD SPACE
FLIGHT CENTER



27 June 2006 Keith Gendreau (GSFC)

ESTC 2006



Prototype Data




GODDARD SPACE FLIGHT CENTER

**Random
Piece of
Quartz
Found
Around
Building 2
at GSFC**

**Clearly see Si
Fluorescence
And ~10 XRD
Signatures.**

**Note that the
arcs are broken
up.**






Prototype Data (Quartz Sample)




GODDARD SPACE
FLIGHT CENTER




2θ

The energy (wavelength) of the events at constant θ are the same. As θ gets smaller, the energy (wavelength) of the events gets higher (shorter) in accordance with Bragg's Law.



Y



X

The combination of physical collecting area and broad wavelength coverage with spectroscopy makes the CXRDF insensitive to sample preparation.